PRODUKTINFORMATION

BESKRIVNING

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* Baron D, Flin C, Porterie J, Despaux J, Vincent P. Hyaluronic acid single intra-articular injection in knee osteoarthritis: a multicenter open prospective study (ART-ONE 75) with placebo post hoc comparison - Current Therapeutic Research 88 (2018) 35-46
<https://doi.org/10.1016/j.curtheres.2018.04.001>

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PRODUCT INFORMATION

DESCRIPTION

ARTHRUM VISC 75 is a sterile viscoelastic, nonpyrogenic and physiological solution of nonreticulated hyaluronic acid of non-animal origin.

This solution comes in a pre-filled single-use syringe. Each box contains a syringe of ARTHRUM VISC 75, instructions and a set of labels to ensure traceability.

COMPOSITION

ARTHRUM VISC 75	3 ml
Sodium hyaluronate (bf)	75 mg
Sodium chloride	27 mg
Borate buffer Ph 7.2 q.s. (Sodium Borate, Boric Acid, Water for Injection)	qsp 3 ml

STERILISATION

The syringes of ARTHRUM VISC 75 are sterilised by moist heat.

INDICATIONS

ARTHRUM VISC 75 is indicated to reduce pain and remobilize the joint by replacing and supplementing the viscoelasticity of the synovial fluid in the arthritic joint.

ARTHRUM VISC 75 treats the joint by viscosupplementation, a procedure to improve the physiological and rheological properties of the synovial fluid in the arthritic knee.

ARTHRUM VISC 75 is indicated to treat the symptoms of osteoarthritis of the knee by intra-articular injection.

ARTHRUM VISC 75 must be injected into the joint by a doctor to relieve pain caused by osteoarthritis of the knee.

CONTRAINDICATIONS

ARTHRUM VISC 75 should not be used for :

- Patients who are hypersensitive to hyaluronic acid.
- Patients with autoimmune diseases, sarcoidosis, multiple severe allergies or cardiovascular disease treated with anticoagulants.
- Patients taking immunosuppressants or interferon.
- Women who are pregnant or breast-feeding.
- Children.

ARTHRUM VISC 75 should not be used for patients with inflammation of the joints (RHEUMATOID ARTHRITIS).

ARTHRUM VISC 75 should not be injected at the same time as an injectable corticosteroid.

PRECAUTIONS FOR USE

ARTHRUM VISC 75 is indicated solely for intraarticular injections.

Do not inject intravascularly.

Do not inject outside the joint cavity or in the synovial membrane or capsule.

INCOMPATIBILITIES

Hyaluronic acid (sodium hyaluronate) solution is incompatible with quaternary ammonium salts such as benzalkonium chloride. ARTHRUM VISC 75 should never come into contact with this type of product or surgical equipment which has been cleaned with this type of product.

ADVERSE REACTIONS

Patients should be informed about potential immediate or delayed adverse reactions in connection with the injection of this product by the practitioner before the procedure.

The various potential adverse reactions are :

- **Immediate adverse reactions related to the injection :**
Pain, redness and swelling at the injection site.
- **Inflammatory effects :**
Pseudogout episode.
Chondrocalcinosis.
Pseudoseptic arthritis.
- **Other side effects :**
The treatment has a low efficacy or little effect in relation to the very advanced stage of the osteoarthritis of the knee according to radiological findings.

New and/or continued side effects should be reported by the patient to the practitioner as soon as possible to remedy the problem using appropriate treatment.

Other unwanted side effects related to the injection must be reported to the distributor and/or manufacturer.

INSTRUCTIONS - DOSAGE

ARTHRUM VISC 75 is to be injected into the joint by an approved specialist practitioner in accordance with applicable local regulations.

The practitioner's technical skills are vital for successful treatment. The device must therefore be used by practitioners who have received specific training in injection techniques. The device must be used in the original packaging provided. No guarantee will be provided should this product be modified or used in any way that is not in accordance with the conditions of use given in these instructions, which may affect its sterility, safety or efficacy.

Practitioners must inform patients about the indications, contraindications, incompatibilities and potential adverse reactions related to the device before starting treatment.

Procedure for use :

- Before injection, the area treated must be carefully disinfected.
 - Ensure that the sterility protector is intact.
 - Take the syringe and remove the protective cover.
 - Holding the main part of the syringe, insert the needle into the Luer lock system.
 - Gently turn the needle clockwise to attach the needle in the Luer lock system.
 - Check to see that the needle is attached.
 - Remove the needle cap.
- Ignoring these instructions can result in the needle becoming detached and/or a leak in the Luer lock.
- The dosage regimen for ARTHRUM VISC 75 is one intra-articular injection.

CLINICAL STUDIES

A published clinical study *, opened and prospective (ART-ONE 75) was carried out on 214 patients with knee osteoarthritis, of which 56% were women, with a Kellgren-Lawrence radiological stage at inclusion varying from I to III

(46% KL III). Their average profile was an age of 62.9 years, a body mass index of 27.2 kg / m² and an anteriority of 4 years for gonarthrosis.

Patients received a single injection of ARTHRUM VISC 75 and were followed up at 30, 60, 120 and 180 days.

As a main criterion, for 199 patients intending to treat, a reduction in pain (WOMAC A index) equal to 28.9 points (base 0-100 points), statistically significant and clinically relevant, was observed at 60 days, by report to the score recorded at the time of inclusion (50.3 points). This improvement continued, reaching 33.6 points at 180 days (gain equal to 67%). In secondary criteria, an improvement was also observed for stiffness (WOMAC B index), function (WOMAC C index), for Quality of life and for disability.

Finally, the percentage of patients responding to treatment according to the OMERACT-OARSI criteria was 86% minimum, from 60 days. The overall tolerance was good, and no serious intra-adverse event was observed in the study, ARTHRUM VISC 75 acting only on the treated joint, without any systemic effect.

Note: In this study, a fraction of the population reviewed at 240 days, showed that the improvement in all the criteria was maintained over time.

* Baron D, Flin C, Porterie J, Despaux J, Vincent P.

Hyaluronic acid single intra-articular injection in knee osteoarthritis: a multicenter open prospective study (ART-ONE 75) with placebo post hoc comparison - Current Therapeutic Research 88 (2018) 35-46
<https://doi.org/10.1016/j.curtheres.2018.04.001>

WARNING

- Check the expiry date on the labelling.
- Do not reuse. Sterility cannot be guaranteed when reused.
- Do not re-sterilise.
- It is recommended that you use sterile Luer-tipped 21G (0.8 x 50 mm) needles and ensure their compatibility with Luer-Lock assemblies before use.
- Used needles and syringes must be disposed of in appropriate containers. Refer to the existing guidelines for disposal.

STORAGE CONDITIONS

- Store between +2°C and +30°C.
- Keep away from light and frost.

ARTHRUM VISC 75 is a CE 0459 medical device, made in France by LCA S.A.

Year of authorisation to affix CE marking : 2009.

Date of review of the product instructions : 2020-05.

For professional use only.

LCA S.A.
9, Allée Prométhée, F-28000 Chartres, France
Phone : +33(0)2 37 33 39 30 Fax : +33(0)2 37 33 39 39
E-mail : lca@lca-pharma.com